

Test Name	Change from LifeLabs Method	Change from Bio Method	BC
a-1Antitrypsin	-10%	none	
Alanine Aminotransferase	none	none	
Albumin	-5%	none	
Alkaline Phos	-7%	none	
Amylase	none	none	
Amylase Urine	none	none	
Angiotensin Converting enzyme	+25%	none	
Apolipoprotein A-1	none	none	
Apolipoprotein B	none	none	
Aspartate Aminotransferase	- 10%	none	
Bicarbonate	- 5%	none	
C Reactive Protein	- 4%	none	
C3 Complement	none	none	
C4 Complement	none	none	
Calcium	none	none	
Calcium Urine	none	none	
Chloride	- 2%	none	
Chloride Urine	none	none	
Cholesterol	none	none	
Citrate	none	none	
Creatine Kinase	none	none	
Creatinine	+ 5%	none	
Creatinine Urine	none	none	
Direct Bilirubin	- 30%	none	
Gamma GT	none	none	
Glucose	none	none	
Glucose Urine	none	none	
Haptoglobin	- 0.1 g/L	none	
HDL Cholesterol	+ 3%	none	
Immunoglobulin A	none	none	
Immunoglobulin G	none	-7%	
Immunoglobulin M	none	-6%	
Inorganic Phosphorus	none	none	
Iron	none	none	
Iron saturation	none	none	
Lactate	none	none	
Lactate Dehydrogenase	- 50%	none	
Lipase	+ 20%	none	
Lipoprotein(a)	none	none	
Lithium	none	none	
Magnesium	none	none	
Potassium	none	none	
Potassium Urine	none	none	
Rhematoid Factor	- 2 IU/mL	-5 IU/mL	

Serum B-HCG	-2 U/L	none
Sodium	none	none
Sodium Urine	none	none
Total Bilirubin	- 20%	none
Total Protein	none	none
Total Protein Urine	none	none
Triglycerides	none	none
UIBC	none	none
Urea	none	none
Urea Urine	none	none
Uric Acid	none	none
Uric Acid Urine	none	none
Urine Albumin	none	none
B-Type Natriuretic Peptide	+ 20 %	none
AM Cortisol	- 20 %	- 44%
PM Cortisol	- 20 %	-44%
Estradiol	- 15%	- 15%
Ferritin	- 15 %	- 12%
Follicle Stimulating Hormone	- 8%	- 5%
Leutenizing Hormone	- 15%	- 4%
Progesterone	- 30%	- 22%
Prolactin	+ 25 %	- 10%
Prostate Specific Antigen	none	none
Thyroid Stimulating Hormone	none	- 24%
Thyroxine, Free	- 20%	- 25%
Triiodothyronine, Free	- 10%	- 20%
Troponin I	Test will be replaced with Troponin T high sensitivity	Test will be replaced with Troponin T high sensitivity
Vitamin B12	none	none