



Patient Instructions

Pinworm

Your health care provider has requested a pinworm specimen.

In order to produce a meaningful laboratory test result, follow these instructions exactly.

1. Holding the cap, remove paddle from the container, do not touch the paddle.
2. Separate the buttocks and press the sticky surface of the paddle against the anal area. The specimen should be obtained early in the morning before arising or bowel movement.
3. Return paddle to container and secure cap tightly.
4. Wash hands thoroughly with soap after you have finished.
5. Print patient information on the label.

Manufacturer Version V302-PIN