












MICROBIOLOGY SWAB AND COLLECTION CONTAINERS

		
<p>Sterile Container (90 mL)</p>	<p>Urine Preservative/ Boric Acid Tube</p>	<p>Pediatric Urine Collector</p>

	
<p>Swabs C&S- Amies Transport Medium Charcoal or Clear Gel</p>	<p>PORT-A-CUL anaerobic transport system</p>

			
<p>Aerobic Blue</p>	<p>Anaerobic Gold</p>	<p>Pediatric Pink</p>	
<p>Bactec Blood Culture Bottles</p>			<p>Fungus Scraping Kit- For skin, hair, or nails</p>

MICROBIOLOGY SWAB AND COLLECTION CONTAINERS

 <p>Molecular testing of CT/GC and TV from:</p> <ul style="list-style-type: none"> • Male Urethral • Female Endocervical & Cervical 	 <p>Molecular testing for the following sources:</p> <ul style="list-style-type: none"> • Vaginal (CT/GC/TV) • Throat (CT & GC) • Rectal (CT & GC)
<p>Aptima Collection Kit</p>	

			
<p>Stool Enteric (Culture & Sensitivity) Cary Blair (Pink liquid)</p>	<p>Stool Ova & Parasites (Manual Microscopy) SAF (Colourless liquid)</p>	<p>Stool Parasites (PCR) Copan Fecal Swab</p>	<p>Pinworm Vials with Paddles</p>

Sources

HOLOGIC, Inc, Aptima® Unisex Swab Specimen Collection Kit Clinician Collection Procedure Guide DS-01696-001 Rev. 003 © 2019
 HOLOGIC, Inc, Aptima® Multitest Swab Specimen Collection Kit Clinician Collection Procedure Guide DS-01552-001 Rev.006 2020
 COPAN, Inc, Copan FecalSwab Collection Kit Collection Procedure Guide HPC021B Rev. 02 Date 2016.11