

Health Care Provider Bulletin

CHANGES TO ANTI-TTG IGA AND ANTI-DGP IGG TESTING

As part of ongoing efforts to improve our quality of service, LifeLabs Ontario is updating the test interpretation for anti-tissue transglutaminase (tTG) IgA antibody and anti-deamidated gliadin peptide (DGP) IgG antibody testing.

This change will go into effect on June 02, 2025.

Table 1. Updates to test interpretation for anti-tTG IgA and anti-DGP IgG

Test	Analyte Concentration	Current Interpretation	Updated Interpretation
Anti-tTG IgA	< 12.0 U/mL	NEGATIVE	NEGATIVE
	12.0 – 14.9 U/mL	BORDERLINE	NEGATIVE
	≥ 15.0 U/mL	POSITIVE	POSITIVE
Anti-DGP IgG	< 12.0 U/mL	NEGATIVE	NEGATIVE
	12.0 – 14.9 U/mL	BORDERLINE	NEGATIVE
	≥ 15.0 U/mL	POSITIVE	POSITIVE

If you have any questions or concerns, please contact the LifeLabs Clinical Biochemist below or LifeLabs Customer Care Centre at 1-877-849-3637.

Dorothy Truong PhD FCACB

Clinical Biochemist

LifeLabs | 100 International Blvd. | Toronto, ON M9W 6J6

T 416-675-4530 Ext. 42208 | C 647-269-7442

E Dorothy.Truong@LifeLabs.com

www.LifeLabs.com