










## MICROBIOLOGY SWAB AND COLLECTION CONTAINERS





		
<p><b>Sterile Container (90 mL)</b></p>	<p><b>Urine Preservative/ Boric Acid Tube</b></p>	<p><b>Pediatric Urine Collector</b></p>
		
<p><b>Swabs C&amp;S- Amies Transport Medium Charcoal or Clear Gel</b></p>		<p><b>E-Swab</b></p>

## MICROBIOLOGY SWAB AND COLLECTION CONTAINERS

			
<p><b>Aerobic Blue</b></p>	<p><b>Anaerobic Lytic</b></p>	<p><b>Pediatric Pink</b></p>	
<p><b>Bactec Blood Culture Bottles</b></p>			<p><b>Fungus Scraping Kit- For skin, hair, or nails</b></p>

	
<p><b>Molecular testing</b> of CT/GC and TV from:</p> <ul style="list-style-type: none"> <li>• <b>Male</b> Urethral</li> <li>• <b>Female</b> Endocervical &amp; Cervical</li> </ul>	<p><b>Molecular testing</b> for the following sources:</p> <ul style="list-style-type: none"> <li>• Vaginal (CT/GC/TV)</li> <li>• Throat (CT &amp; GC)</li> <li>• Rectal (CT &amp; GC)</li> </ul>
<p><b>Aptima Collection Kit</b></p>	

## MICROBIOLOGY SWAB AND COLLECTION CONTAINERS

			
<p><b>Stool Enteric (Culture &amp; Sensitivity)</b> <b>Cary Blair (Pink liquid)</b></p>	<p><b>Stool Ova &amp; Parasites (Manual Microscopy)</b> <b>SAF (Colourless liquid)</b></p>	<p><b>Stool Parasites (PCR)</b> <b>Copan Fecal Swab</b></p>	<p><b>Pinworm Vials with Paddles</b></p>

**Sources**

HOLOGIC, Inc, Aptima® Unisex Swab Specimen Collection Kit Clinician Collection Procedure Guide DS-01696-001 Rev. 003 © 2019  
 HOLOGIC, Inc, Aptima® Multitest Swab Specimen Collection Kit Clinician Collection Procedure Guide DS-01552-001 Rev.006 2020  
 COPAN, Inc, Copan FecalSwab Collection Kit Collection Procedure Guide HPC021B Rev. 02 Date 2016.11